

CareSens[®] N Eco

Blood Glucose Monitoring System



Display with Intuitive Symbols

- Customizable range indicator
- Small sample size with rapid test result



Specifications

- 0.5 μ L / 5 seconds
- 20-600 mg/dL (1.1-33.3 mmol/L)
- Hematocrit: 15-65%
- 1,000 test result memory
- 5-45°C (41-113°F)

Features

- No coding required
- Strip ejection button
- Normal range, Hypoglycemia/Hyperglycemia indicator
- Strip expiration date indicator
- Beep sound on/off function
- 4 alarms (PP2* and 3 time-set alarms)
*2 hours after meal
- 1, 7, 14, 30, 90-day test averages
- Pre-meal, Post-meal, Fasting, No-mark flagging
- Data port

Easy Diabetes Care: Less work and less worry with CareSens® N Eco

One-step Operation

Just insert the strip and you are ready to test in one step



① Insert test strip






② Apply blood sample



③ Result in 5 seconds

Range Indicator with Intuitive Symbols

- Identify test results with pre-set ranges through symbol display
- Desired range for each indicator can be customized



LOW

Hypoglycemia Indicator



IN RANGE

Normal range Indicator



HIGH

Hyperglycemia Indicator

Improved Measurement Performance

- Advanced accuracy
 - 99.0% of results within ± 15 mg/dL (0.83 mmol/L) or 15% of the reference results, meeting ISO 15197:2013 (EN ISO 15197:2015) accuracy criteria

Within ± 5 mg/dL or $\pm 5\%$	Within ± 10 mg/dL or $\pm 10\%$	Within ± 15 mg/dL or $\pm 15\%$
389/600 (64.8%)	552/600 (92.0%)	389/600 (99.0%)

